

EBP05417 - IMMUNIZING PEPTIDE

SKU: EBP05417



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_005301.2
Immunogen	Peptide with sequence ELDLSPPHLSSSPE-C, from the N Terminus of the protein sequence according to NP_005301.2.
Product Comments	The peptide was used in the production of EB05417 - Goat Anti-GRB7 (N Terminus) Antibody. It is currently untested in blocking applications.
Peptide Sequence	ELDLSPPHLSSSPE-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)