

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05419-B - IMMUNIZING PEPTIDE

SKU: EBP05419-B

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 116573.1; NP 004838.1; NP 001614.3 ID

Peptide with sequence C-TGPPAKKAISELP., from the C Terminus of the protein sequence according to Immunogen

NP_116573.1; NP_004838.1; NP_001614.3.

The peptide was used in the production of EB05419-B - Goat Anti-AIF1 / IBA1 (isoform 1 and 3), biotinylated **Product**

Comments Antibody. It is currently untested in blocking applications.

Peptide C-TGPPAKKAISELP. Sequence

Shipping

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade