

# EBP05419-B - IMMUNIZING PEPTIDE

**SKU:** EBP05419-B



---

## SPECIFICATIONS

|                              |   |
|------------------------------|---|
| <b>Formulation</b>           | 100µg of Lyophilised peptide, formulation De-ionized (DI) water.  |
| <b>Unit Size</b>             | 100 µg  |
| <b>Storage</b>               | Shipped at ambient temperature, store at -20°C.   |
| <b>Instructions</b>          |   |
| <b>Accession ID</b>          | NP_116573.1; NP_004838.1; NP_001614.3   |
| <b>Immunogen</b>             | Peptide with sequence C-TGPPAKKAISELP., from the C Terminus of the protein sequence according to NP_116573.1; NP_004838.1; NP_001614.3.   |
| <b>Product</b>               | The peptide was used in the production of EB05419-B - Goat Anti-AIF1 / IBA1 (isoform 1 and 3), biotinylated Antibody. It is currently untested in blocking applications.  |
| <b>Comments</b>              |   |
| <b>Peptide Sequence</b>      | C-TGPPAKKAISELP.  |
| <b>Shipping Instructions</b> | Ambient   |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a> |

## DOCUMENTS

- [Data Sheet](#)