

EBP05419-H - IMMUNIZING PEPTIDE

SKU: EBP05419-H



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_116573.1; NP_001614.3
Immunogen	Peptide with sequence C-TGPPAKKAISELP.., from the C Terminus of the protein sequence according to NP_116573.1; NP_001614.3.
Product	The peptide was used in the production of EB05419-H - Goat Anti-AIF1/IBA1 isoform 1 and 3 (Human reactivity only) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-TGPPAKKAISELP..
Shipping	
Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png