

EBP05419-R - IMMUNIZING PEPTIDE

SKU: EBP05419-R



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_116573.1; NP_001614.3

Immunogen Peptide with sequence C-TGPPAKKAISELP.,, from the C Terminus of the protein sequence according to NP_116573.1; NP_001614.3.

Product Comments The peptide was used in the production of EB05419-R - Goat Anti-AIF1/IBA1 isoform 1 and 3 (Rat reactivity only) Antibody. It is currently untested in blocking applications.

Peptide Sequence C-TGPPAKKAISELP.,

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png