



EBP05445 - IMMUNIZING PEPTIDE

SKU: EBP05445



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_065122.1; NP_001161049.1; NP_001161048.1
Immunogen	Peptide with sequence C-NKDHLRVNKGKDI, from the C Terminus of the protein sequence according to NP_065122.1; NP_001161049.1; NP_001161048.1.
Product Comments	The peptide was used in the production of EB05445 - Goat Anti-TRP7 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-NKDHLRVNKGKDI
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png