



EBP05447 - IMMUNIZING PEPTIDE

SKU: EBP05447



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_057399.1
Immunogen	Peptide with sequence EGTVFCLDPLDSRC, from the C Terminus of the protein sequence according to NP_057399.1.
Product	The peptide was used in the production of EB05447 - Goat Anti-GULP1 / CED6 Antibody. It is currently untested
Comments	in blocking applications.
Peptide Sequence	EGTVFCLDPLDSRC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)