

EBP05465 - IMMUNIZING PEPTIDE

SKU: EBP05465



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_954870.2; NP_954868.1; NP_954869.1
Immunogen	Peptide with sequence C-EKMTSSDVVAGSD, from the C Terminus of the protein sequence according to NP_954870.2; NP_954868.1; NP_954869.1.
Product Comments	The peptide was used in the production of EBP05465 - Goat Anti-TPTE / PTEN2 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-EKMTSSDVVAGSD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png