

EBP05476 - IMMUNIZING PEPTIDE

SKU: EBP05476



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_443123.1; NP_065126.2
Immunogen	Peptide with sequence C-ALYNIISTWPHFKH, from the C Terminus of the protein sequence according to NP_443123.1; NP_065126.2.
Product	The peptide was used in the production of EB05476 - Goat Anti-PGLYRP3 / PGLYRP4 Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	C-ALYNIISTWPHFKH
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png