

EBP05487 - IMMUNIZING PEPTIDE

SKU: EBP05487



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001015049.1; NP_004864.1; NP_001015048.1
Immunogen	Peptide with sequence DMGNQHPSISRLQ-C, from the N Terminus of the protein sequence according to NP_001015049.1; NP_004864.1; NP_001015048.1.
Product Comments	The peptide was used in the production of EB05487 - Goat Anti-BAG5 (N Terminus) Antibody. It is currently untested in blocking applications.
Peptide Sequence	DMGNQHPSISRLQ-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)