

EBP05488 - IMMUNIZING PEPTIDE

SKU: EBP05488



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_689452.1; NP_689451.1; NP_004368.1; NP_001138609.1; NP_001138607.1; NP_001138606.1; NP_001138608.1
Immunogen	Peptide with sequence C-DIADLFQVPKAYS, from the C Terminus of the protein sequence according to NP_689452.1; NP_689451.1; NP_004368.1; NP_001138609.1; NP_001138607.1; NP_001138606.1; NP_001138608.1.
Product Comments	The peptide was used in the production of EB05488 - Goat Anti-CPT1B (isoform 1) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-DIADLFQVPKAYS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)