

EBP05502 - IMMUNIZING PEPTIDE

SKU: EBP05502



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_004865.1

Immunogen Peptide with sequence SALRRSGYGPSDGP-C, from the N Terminus of the protein sequence according to NP_004865.1.

Product Comments The peptide was used in the production of EB05502 - Goat Anti-BAG4 / SODD Antibody. It is currently untested in blocking applications.

Peptide Sequence SALRRSGYGPSDGP-C

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png