

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05505 - IMMUNIZING PEPTIDE

SKU: EBP05505

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_057576.2 ID

Immunogen Peptide with sequence C-SLTLHDPNIFGPQ, from the C Terminus of the protein sequence according to NP_057576.2.

The peptide was used in the production of EB05505 - Goat Anti-RANGNRF / MOG1 Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-SLTLHDPNIFGPQ Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade