

EBP05511 - IMMUNIZING PEPTIDE

SKU: EBP05511



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_938195.1; NP_001413.1, NP_001230009.1, NP_001230010.1
Immunogen	Peptide with sequence C-GKNAHWQEDKL, from the C Terminus of the protein sequence according to NP_938195.1; NP_001413.1, NP_001230009.1, NP_001230010.1.
Product Comments	The peptide was used in the production of EB05511 - Goat Anti-ELF5 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-GKNAHWQEDKL
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiling_medium.png

DOCUMENTS

- [Data Sheet](#)