

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05511 - IMMUNIZING PEPTIDE

SKU: EBP05511

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_938195.1; NP_001413.1, NP_001230009.1, NP_001230010.1 ID

The peptide was used in the production of EB05511 - Goat Anti-ELF5 Antibody. It is currently untested in blocking Product

Comments applications.

Peptide

C-GKNAHGWQEDKL Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade