

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05513 - IMMUNIZING PEPTIDE

SKU: EBP05513

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 003457.1; NP 039245.1; NP 039246.1; NP 039247.1; NP 054698.1 ID

Peptide with sequence PHNSIRSGHGGLN-C, from the N Terminus of the protein sequence according to Immunogen

NP 003457.1; NP 039245.1; NP 039246.1; NP 039247.1; NP 054698.1.

The peptide was used in the production of EB05513 - Goat Anti-PAX8 (N Terminus) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

PHNSIRSGHGGLN-C Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png