

EBP05530 - IMMUNIZING PEPTIDE

SKU: EBP05530



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001035197.1; NP_055196.2
Immunogen	Peptide with sequence RQSLQQSMPPCEPP, from the C Terminus of the protein sequence according to NP_001035197.1; NP_055196.2.
Product	The peptide was used in the production of EB05530 - Goat Anti-MLH3 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	RQSLQQSMPPCEPP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)