

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05539 - IMMUNIZING PEPTIDE

SKU: EBP05539

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_000525.1; NP_001121615.1; NP_001121616.1 ID

The peptide was used in the production of EB05539 - Goat Anti-PMS1 Antibody. It is currently untested in blocking Product

Comments applications.

Peptide

C-RPFFHHLTYLPETT Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade