

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05542 - IMMUNIZING PEPTIDE

SKU: EBP05542

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 001168 ID

Peptide with sequence C-EAMEWLVETLKSRQ, from the C Terminus of the protein sequence according to **Immunogen**

NP_001168.

The peptide was used in the production of EB05542 - Goat Anti-ARL1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide C-EAMEWLVETLKSRQ Sequence

Shipping

Instructions Ambient

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$ Grade