

EBP05547 - IMMUNIZING PEPTIDE

SKU: EBP05547



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001652.2
Immunogen	Peptide with sequence GNHITEMAPTASS-C, from the N Terminus of the protein sequence according to NP_001652.2.
Product Comments	The peptide was used in the production of EB05547 - Goat Anti-ARL4D / ARL4L Antibody. It is currently untested in blocking applications.
Peptide Sequence	GNHITEMAPTASS-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png