

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05549 - IMMUNIZING PEPTIDE

SKU: EBP05549

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 001032241.1; NP 005729.1; P 997625.1 ID

Peptide with sequence GNGLSDQTSILSN-C, from the N Terminus of the protein sequence according to Immunogen

NP 001032241.1; NP 005729.1; P 997625.1.

The peptide was used in the production of EB05549 - Goat Anti-ARL4A Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide GNGLSDQTSILSN-C

Sequence Shipping

Product

Instructions Ambient

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade