

EBP05559 - IMMUNIZING PEPTIDE

SKU: EBP05559



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_006568; NP_150274

Immunogen Peptide with sequence DIWSQLQSAHLKC, from the C Terminus of the protein sequence according to NP_006568; NP_150274.

Product Comments The peptide was used in the production of EBP05559 - Goat Anti-B3GNT1 Antibody. It is currently untested in blocking applications.

Peptide Sequence DIWSQLQSAHLKC

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png