

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP05569 - IMMUNIZING PEPTIDE**

**SKU:** EBP05569

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

 $NP\_008950.1;\ NP\_861515.1;\ NP\_861516.1;\ NP\_861517.1;\ NP\_861518.1;\ NP\_001268670.1$ ID

Peptide with sequence C-QETYSKQVTSQEP, from the C Terminus of the protein sequence according to NP\_008950.1; Immunogen NP\_861515.1; NP\_861516.1; NP\_861517.1; NP\_861518.1; NP\_001268670.1.

The peptide was used in the production of EB05569 - Goat Anti-UBE2C / UBCH10 Antibody. It is currently untested in **Product** 

Comments blocking applications.

**Peptide** 

C-QETYSKQVTSQEP Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade