

EBP05584 - IMMUNIZING PEPTIDE

SKU: EBP05584



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_057168.1
Immunogen	Peptide with sequence GSQAKDKERCVIS, from the C Terminus of the protein sequence according to NP_057168.1.
Product	The peptide was used in the production of EBP05584 - Goat Anti-RASD1 / AGS1 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	GSQAKDKERCVIS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png