

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP05587 - IMMUNIZING PEPTIDE**

**SKU:** EBP05587

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

100 μg **Unit Size** 

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 060469.4 ID

Peptide with sequence C-EEVPLGQSIPRRRK, from the C Terminus of the protein sequence according to **Immunogen** 

NP\_060469.4.

The peptide was used in the production of EB05587 - Goat Anti-RNF31 / ZIBRA (C terminus) Antibody. It is **Product** 

**Comments** currently untested in blocking applications.

**Peptide** 

C-EEVPLGQSIPRRRK Sequence

**Shipping** 

Instructions Ambient

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png