

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05619 - IMMUNIZING PEPTIDE

SKU: EBP05619

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 001021.1 ID

Peptide with sequence PLAKDLLHPSPEE-C, from the N Terminus of the protein sequence according to **Immunogen**

NP 001021.1.

The peptide was used in the production of EB05619 - Goat Anti-MPS1 / RPS27 Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

PLAKDLLHPSPEE-C Sequence

Shipping Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade