



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05631 - IMMUNIZING PEPTIDE

SKU: EBP05631

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_006398.1 ID

Immunogen Peptide with sequence C-NRLTDYISKVKE, from the C Terminus of the protein sequence according to NP_006398.1. **Product** The peptide was used in the production of EB05631 - Goat Anti-GTRAP3-18 / JWA Antibody. It is currently untested

Comments in blocking applications.

Peptide

C-NRLTDYISKVKE Sequence

Shipping

Ambient Instructions

Grade

Product $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$