

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05674-B - IMMUNIZING PEPTIDE

SKU: EBP05674-B

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 000007.1; NP 001120800.1; NP 001272971.1; NP 001272972.1; NP 001272973.1 ID

Peptide with sequence C-RLIVAREHIDKYKN., from the C Terminus of the protein sequence according to Immunogen

NP 000007.1; NP 001120800.1; NP 001272971.1; NP 001272972.1; NP 001272973.1.

The peptide was used in the production of EB05674-B - Goat Anti-ACADM, Biotinylated Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-RLIVAREHIDKYKN. Sequence

Shipping

Instructions Ambient

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png