

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP05685 - IMMUNIZING PEPTIDE**

**SKU:** EBP05685

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

100 μg **Unit Size** 

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 829883.1; NP 829884.1 ID

Peptide with sequence YGSARSVGKVEPS-C, from the N Terminus of the protein sequence according to Immunogen

NP 829883.1; NP 829884.1.

The peptide was used in the production of EB05685 - Goat Anti-ELKS / RAB6IP2 Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

YGSARSVGKVEPS-C Sequence

Shipping

Instructions Ambient

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade