

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## EBP05687 - IMMUNIZING PEPTIDE

**SKU:** EBP05687

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession NP\_002868.1; NP\_598194.1; NP\_598193.2; NP\_001308738.1; NP\_001308741.1; NP\_001308747.1;

NP 001308748.1;; NP 001308750.1 ID

Peptide with sequence GSKKLKRVGLSQE-C, from the N Terminus of the protein sequence according to

Immunogen NP\_002868.1; NP\_598194.1; NP\_598193.2; NP\_001308738.1; NP\_001308741.1; NP\_001308747.1;

NP\_001308748.1;; NP\_001308750.1.

The peptide was used in the production of EB05687 - Goat Anti-RAD51L1 Antibody. It is currently untested in **Product** 

Comments blocking applications.

**Peptide** 

GSKKLKRVGLSQE-C Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade