

EBP05687 - IMMUNIZING PEPTIDE

SKU: EBP05687



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002868.1; NP_598194.1; NP_598193.2; NP_001308738.1; NP_001308741.1; NP_001308747.1; NP_001308748.1;; NP_001308750.1
Immunogen	Peptide with sequence GSKKLKRVGLSQE-C, from the N Terminus of the protein sequence according to NP_002868.1; NP_598194.1; NP_598193.2; NP_001308738.1; NP_001308741.1; NP_001308747.1; NP_001308748.1;; NP_001308750.1.
Product Comments	The peptide was used in the production of EB05687 - Goat Anti-RAD51L1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	GSKKLKRVGLSQE-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)