

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP05691 - IMMUNIZING PEPTIDE**

**SKU:** EBP05691

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_001124334.1; NP\_006470 ID

Immunogen Peptide with sequence C-ARVKPLHPNATST, from the C Terminus of the protein sequence according to NP\_001124334.1; NP\_006470.

**Product** The peptide was used in the production of EB05691 - Goat Anti-PIR51 / RAD51AP1 Antibody. It is currently untested

**Comments** in blocking applications.

**Peptide** 

C-ARVKPLHPNATST Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite\_plus\_medium.png$ Grade