

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05692 - IMMUNIZING PEPTIDE

SKU: EBP05692

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_002869.3; NP_598332.1; NP_001136043.1 ID

The peptide was used in the production of EB05692 - Goat Anti-RAD51D / RAD51L3 Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

GVLRVGLCPGLTEE Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade