

EBP05692 - IMMUNIZING PEPTIDE

SKU: EBP05692



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002869.3; NP_598332.1; NP_001136043.1
Immunogen	Peptide with sequence GVLRVGLCPGLTEE, from the N Terminus of the protein sequence according to NP_002869.3; NP_598332.1; NP_001136043.1.
Product Comments	The peptide was used in the production of EBP05692 - Goat Anti-RAD51D / RAD51L3 Antibody. It is currently untested in blocking applications.
Peptide Sequence	GVLRVGLCPGLTEE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)