

EBP05696 - IMMUNIZING PEPTIDE

SKU: EBP05696



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_036300.2
Immunogen	Peptide with sequence C-GARVTNSSVWVEP, from the C Terminus of the protein sequence according to NP_036300.2.
Product Comments	The peptide was used in the production of EB05696 - Goat Anti-FBX2 / FBXO2 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-GARVTNSSVWVEP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png