

EBP05714 - IMMUNIZING PEPTIDE

SKU: EBP05714



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_061169.2
Immunogen	Peptide with sequence C-KVHMDNIHGVADS, from the C Terminus of the protein sequence according to NP_061169.2.
Product Comments	The peptide was used in the production of EB05714 - Goat Anti-PRDM5 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KVHMDNIHGVADS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png