

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05717 - IMMUNIZING PEPTIDE

SKU: EBP05717

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_064613.2; NP_955469.1; NP_955470.1; NP_955471.1; NP_001354824.1; NP_001354826.1 ID

Peptide with sequence C-NGNGSSEVHITKP, from the C Terminus of the protein sequence according to

 $NP_064613.2;\ NP_955469.1;\ NP_955470.1;\ NP_955471.1;\ NP_001354824.1;\ NP_001354826.1.$

The peptide was used in the production of EB05717 - Goat Anti-PRDM10 Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

C-NGNGSSEVHITKP Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade