

EBP05722 - IMMUNIZING PEPTIDE

SKU: EBP05722



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_006700; NP_079532

Immunogen Peptide with sequence C-LPELGSLPPVNT, from the C Terminus of the protein sequence according to NP_006700; NP_079532.

Product Comments The peptide was used in the production of EB05722 - Goat Anti-BAT8 / G9A Antibody. It is currently untested in blocking applications.

Peptide Sequence C-LPELGSLPPVNT

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png