

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05738 - IMMUNIZING PEPTIDE

SKU: EBP05738

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

 $NP_001073289.1; \ NP_004886.3; \ NP_899632.1; \ NP_001120933.1; \ NP_001120934.1$ ID

Peptide with sequence C-EKPELTVVFEPSW, from the C Terminus of the protein sequence according to

Immunogen NP_001073289.1; NP_004886.3; NP_899632.1; NP_001120933.1; NP_001120934.1.

The peptide was used in the production of EB05738 - Goat Anti-ClAS1 / cryopyrin (C Terminus) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-EKPELTVVFEPSW Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade