

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05750 - IMMUNIZING PEPTIDE

SKU: EBP05750

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

ID

NP 005330.1; NP 001104582.1; NP 001104583.1

Peptide with sequence C-SWDVETATELLLSN, from the C Terminus of the protein sequence according to Immunogen

NP 005330.1; NP 001104582.1; NP 001104583.1.

The peptide was used in the production of EB05750 - Goat Anti-HIP2 Antibody. It is currently untested in **Product**

Comments blocking applications. **Peptide**

C-SWDVETATELLLSN Sequence

Shipping Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade