

EBP05761 - IMMUNIZING PEPTIDE

SKU: EBP05761



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_071349.3

Immunogen Peptide with sequence ALLEAGMPECTEDK, from the C Terminus of the protein sequence according to NP_071349.3.

Product Comments The peptide was used in the production of EB05761 - Goat Anti-E2-230K / UBE2O Antibody. It is currently untested in blocking applications.

Peptide Sequence ALLEAGMPECTEDK

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)