

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP05765 - IMMUNIZING PEPTIDE**

**SKU:** EBP05765

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_054895.1; NP\_001297255.1 ID

 $\label{eq:lower_$ 

The peptide was used in the production of EB05765 - Goat Anti-HSPC150 / UBE2T Antibody. It is currently untested **Product** 

Comments in blocking applications.

**Peptide** 

C-QLVGIEKKFHPDV Sequence

**Shipping** 

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite\_plus\_medium.png$ Grade