

EBP05774 - IMMUNIZING PEPTIDE

SKU: EBP05774



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_055107.3; NP_065700.1; NP_001290423.1
Immunogen	Peptide with sequence ELSAVGERVFAAE-C, from the N Terminus of the protein sequence according to NP_055107.3; NP_065700.1; NP_001290423.1.
Product Comments	The peptide was used in the production of EB05774 - Goat Anti-CBX6 and CBX8 Antibody. It is currently untested in blocking applications.
Peptide Sequence	ELSAVGERVFAAE-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)