

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP05798 - IMMUNIZING PEPTIDE**

**SKU:** EBP05798

## **SPECIFICATIONS**

Formulation  $100\mu g$  of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μα

Shipped at ambient temperature, store at -20°C.

Accession

NP\_006615.2; NP\_997644.2; NP\_001189393.1; NP\_001189394.1; NP\_001189395.1; NP\_001189396.1; NP\_001189397.1 ID

Peptide with sequence SRVHGMHPKETT-C, from the N Terminus of the protein sequence according to

Immunogen NP 006615.2; NP 997644.2; NP 001189393.1; NP 001189394.1; NP 001189395.1; NP 001189396.1;

NP 001189397.1.

The peptide was used in the production of EB05798 - Goat Anti-BS69 / ZMYND11 Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

SRVHGMHPKETT-C Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade