

EBP05802 - IMMUNIZING PEPTIDE

SKU: EBP05802



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_037508.1

Immunogen Peptide with sequence LSKGLKRKREEEE-C, from the N Terminus of the protein sequence according to NP_037508.1.

Product Comments The peptide was used in the production of EB05802 - Goat Anti-SEI-1 / TRIP-Br1 Antibody. It is currently untested in blocking applications.

Peptide Sequence LSKGLKRKREEEE-C

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png