

# EBP05807 - IMMUNIZING PEPTIDE

**SKU:** EBP05807



---

## SPECIFICATIONS

|                              |   |
|------------------------------|---|
| <b>Formulation</b>           | 100µg of Lyophilised peptide, formulation De-ionized (DI) water.  |
| <b>Unit Size</b>             | 100 µg  |
| <b>Storage</b>               | Shipped at ambient temperature, store at -20°C.   |
| <b>Instructions</b>          |   |
| <b>Accession ID</b>          | NP_056989.2; NP_003843.3  |
| <b>Immunogen</b>             | Peptide with sequence C-KRLKSIEERQLLK, from the C Terminus of the protein sequence according to NP_056989.2; NP_003843.3.   |
| <b>Product Comments</b>      | The peptide was used in the production of EB05807 - Goat Anti-TIF1A / TRIM24 Antibody. It is currently untested in blocking applications.   |
| <b>Peptide Sequence</b>      | C-KRLKSIEERQLLK   |
| <b>Shipping Instructions</b> | Ambient   |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a> |

## DOCUMENTS

- [Data Sheet](#)