

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05818 - IMMUNIZING PEPTIDE

SKU: EBP05818

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_006559.1 ID

Immunogen Peptide with sequence CSQKEQDKDKPKTL, from the C Terminus of the protein sequence according to NP_006559.1.

The peptide was used in the production of EB05818 - Goat Anti-CGR19 / CGRRF1 Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

CSQKEQDKDKPKTL Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade