

EBP05819 - IMMUNIZING PEPTIDE

SKU: EBP05819



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001154816.1; NP_001154817.1; NP_001154818.1; NP_060693.2; NP_001154819.1
Immunogen	Peptide with sequence KFNHICEQTRFKN, from the C Terminus of the protein sequence according to NP_001154816.1; NP_001154817.1; NP_001154818.1; NP_060693.2; NP_001154819.1.
Product Comments	The peptide was used in the production of EB05819 - Goat Anti-CHFR Antibody. It is currently untested in blocking applications.
Peptide Sequence	KFNHICEQTRFKN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)