

EBP05825 - IMMUNIZING PEPTIDE

SKU: EBP05825



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_004632.1; NP_001009569.1

Immunogen Peptide with sequence C-HFSPSMEIMQVRK, from the C Terminus of the protein sequence according to NP_004632.1; NP_001009569.1.

Product Comments The peptide was used in the production of EBP05825 - Goat Anti-AF10 / MLLT10 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-HFSPSMEIMQVRK

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png