

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **EBP05831 - IMMUNIZING PEPTIDE**

**SKU:** EBP05831

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

ID

NP\_006596.2; NP\_001091997.2; NP\_001153017.1; NP\_001153018.1; NP\_001092005.1; NP\_006595.1

Peptide with sequence C-TTDAPVRPGEAK, from the C Terminus of the protein sequence according to NP\_006596.2;

 $NP\_001091997.2;\ NP\_001153017.1;\ NP\_001153018.1;\ NP\_001092005.1;\ NP\_006595.1.$ 

The peptide was used in the production of EB05831 - Goat Anti-RFPL2 and RFPL3 Antibody. It is currently untested Product

**Comments** in blocking applications.

**Peptide** 

C-TTDAPVRPGEAK Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade